Innovint Aircraft Interior GmbH

entwicklung von innenausstattung innovation of interior Innovation d'équipements commerciaux



ATA CHAPTER





- Airbus
- Boeing
- Embraer
- → ATR
- ... Maintenance MRO's

- **Emergency Medical Equipment**
- → FAR-121 Subpart X
- → EASA AIR-OPS
- → NAA Requirements
- Customer Requirements +++ A.O.G.Service +++ Consulting

First Aid **Kits**

Emergency Medical Kits

Replenishment **Packages**

Customized

Medical

Equipment

Doctor's

Kits

Spare Parts Support





Overview Emergency Medical Equipment

- → BASIC First Aid Kit
- Single Aisle First Aid Kit
- Long Range First Aid Kit
- → EASA AIR-OPS Emergency Medical Kit
- → FAA Emergency Medical Kit
- Universal Precaution Kit (Standard or Long Range version)
- Attendant Kit
- In case of other needed requirements Innovint can provide a compliance analysis.

The medical equipment is defined by the FAR, EASA and ICAO requirements.

The requirements prescribe the minimum mandatory standard. It is always the airline decision to select advanced configurations regarding contents in order to provide cabin crew and voluntary physicians with a wider range of equipment, thus saving many lives in an emergency medical situation. Innovint kits are based on a modular system and comprise mandatory basic kits as well as advanced and customized configurations.



PNR 6240-series • First Aid Kit stored in bracket

EASA AIR-OPS REQUIREMENTS

First Aid Kit EAS	DA AIN-OF S
1-99 pax 1 ea. 100-199 pax 2 ea. 200-299 pax 3 ea. 300 up pax 4 ea.	

Emergency Med. Kit EASA AIR-OPS

Required for more than 30 pax. if any point on the planned route is more than 60 minutes flying time from an airport at which qualified medical assistance is assured.



PNR 6240-series • Emergency Medical Kit

FAR-121 Subpart X

First Aid K	<i>Tit</i>		FAR-121
50 51-150 151-250 250 up	pax pax pax pax	1 ea. 2 ea. 3 ea. 4 ea.	
Emergenc	y Med. k	(it	FAR-121

One approved emergency medical kit in aircraft for which a flight attendant is required.

ed. 06/15 subject to alterations





BASIC First Aid Kit

- Airbus cross fleet approval
- Boeing basic approval
- → ETSO Approval under technical specification TS 40.400/2TS
- → In compliance with EASA AIR-OPS
- → In compliance with FAR-121, subpart X
- → In compliance with ICAO, part I

Anticipated shelf life is 18 months, which may be reduced down to 12 months due to the availability of medicine production lots. Latches can be sealed by tie-wraps, plombs or lockers. The first aid kit PNR 6240-37 has been strictly assembled to the statutory requirements of EASA AIR-OPS. In addition the kits cover FAR 121, subpart X and ICAO, part I as well. It represents the classical standard equipment. The contents of first aid kit PNR 6240-48 is reduced for FAR-121 requirements only. The kits are manufactured of aluminum sheets (coated in grey) and aluminum sections.

The bracket PNR 6240-01-001 allows installation of the kit at any location. The bracket has been tested in our approved testing facilities for use in all directions (9g).



DIMENSIONS

First Aid Kit

-	PNR 6240-37 /-48


First Aid Kit




Length:	270 mm	10.63 in
Width:	300 mm	11.81 in
Height:	115 mm	4.53 in
Weight max: 6240-37/-67	3 ka	6.61 lbs

2.5 kg

Weight max: 6240-48 (incl. contents)





DIMENSIONS				
Bracket	PNR (6240-37-011		
Length: Width:	278 mm 130 mm	11.00 in 5.12 in		
Height:	280 mm	11.02 in		
Weight:	1 kg	2.20 lbs		

PNR 6240-37 • Bracket

ed. 04/15 subject to alterations

5.50 lbs





Single Aisle First Aid Kit

- Airbus cross fleet approval
- Boeing basic approval
- → ETSO Approval under technical specification TS 40.400/2TS
- → In compliance with EASA AIR-OPS
- → In compliance with FAR-121, subpart X
- In compliance with ICAO, part I
- Content exceeds named regulations


Anticipated shelf life is 18 months, which may be reduced down to 12 months due to the availability of medicine production lots. Latches can be sealed by tie-wraps, plombs or lockers.



PNR 6240-35 • First Aid Kit standard body

Thirtoz-45 55 * Finst-And Nit Standard Body

PNR 6240–45 • White powder – coated



These First Aid Kits have been assembled to the statutory requirements of EASA AIR-OPS. In addition the kits cover FAR 121, subpart X and ICAO, part I as well. They are especially designed to be stored in the overhead stowage compartment, taking as little space as possible. PNR 6240-35/-65 is manufactured of aluminum sheets (coated in grey) and aluminum sections.


PNR 6240-45/-46 is made of a more light weight box (coated in white).

The bracket PNR 6240-05-001 allows installation of the kit at any location. It has been tested in our approved testing facilities for use in all directions (9g).

an conone (cg).			
DIMENSIONS			
First Aid Kit		PNR	6240-35 /-45
First Aid Kit	C	PNR	6240-65 /-46
Length:	5	525 mm	20.67 in
Width:	2	205 mm	8.10 in
Height:	1	15 mm	4.53 in
Weight max: 6240-35	5/-65	3.9 kg	8.59 lbs
Weight max: 6240-45	5/-46	3.4 kg	7.49 lbs



DIMENCIONO



Airbus standard location in the overhead stowage compartment

DIMENSIONS				
Bracket	P	NR 6240-05-001		
Length:	506 mm	20.00 in		
Width:	126 mm	5.00 in		
Height:	210 mm	8.27 in		
Weight:	1 kg	2.20 lbs		
		ed. 04/15 subject to alterations		





Long Range First Aid Kit

- Airbus cross fleet approval
- → Boeing basic approval
- → ETSO Approval under technical specification TS 40.400/2TS
- → In compliance with EASA AIR-OPS
- → In compliance with FAR-121, subpart X
- → In compliance with ICAO, part I
- Content exceeds named regulations

Anticipated shelf life is 18 months, which may be reduced down to 12 months due to the availability of medicine production lots. Latches can be sealed by tie-wraps, plombs or lockers. PNR 6240-31 long range first aid kit gives high degree of flexibility for contents customization. Half of the size is required to comply with the EASA AIR-OPS, FAR-121 requirements and ICAO. The other half can be equipped to the airliner needs. The -xx PNR will change accordingly (PNR 6240-31 shown). The kit is manufactured of aluminum sheets (coated in grey) and aluminum sections.

PNR 6240-41 is made of a more light weight box (coated in white).

The bracket PNR 6240-01-001 allows installation of the kits at any location. The bracket has been tested in our approved testing facilities for use in all directions (9g).



PNR 6240-31 • First Aid Kit long range

DIMENSIONS		
First Aid Kit		PNR 6240-31/-41
First Aid Kit	C	PNR 6240-xx/-42
Length:	540 mr	m 21.26 in
Width:	300 mr	n 11.81 in
Height:	115 mr	m 4.53 in
Weight max: 6240-31	4.8 k	g 10.58 lbs
Weight max: 6240-41/-4 (incl. contents)	12 4.2 k	g 9.25 lbs



PNR 6240-41 • White powder - coated

DIMENSIONS				
Bracket	PN	R 6240-01-001		
Length: Width: Height:	551 mm 130 mm 285 mm	21.69 in 5.12 in 11.22 in		
Weight:	1,2 kg	2.64 lbs		

ed. 04/15 subject to alterations



ATA CHAPTER

Soft First Aid Kit

- → In compliance with EASA AIR-OPS
- → In compliance with FAR-121, subpart X
- → In compliance with ICAO, part I

Anticipated shelf life is 18 months, which may be reduced down to 12 months due to the availability of medicine production lots.

Zipper can be sealed by tie-wraps, plombs or lockers.



PNR 6240-87S • First Aid Kit, closed

The first aid kit PNR 6240-87S has been strictly assembled to the statutory requirements of EASA AIR-OPS. In addition the kit covers FAR 121, subpart X and ICAO, part I. The red soft bag is made from robust fire retardant material. The first aid kit is especially designed to be stored in AIRBUS A320 headrests of cabin attendant seats, hatrack lifeline spacers and several other unused cabin stowage spaces. It is light weight, easy to stow and has a transparent content paper card pocket visualizing relevant information like serial number and expiry date.

DIMENSIONS				
First Aid Kit		PNR 6240-87S		
Length:	280 mm	11.02 in		
Width:	280 mm	11.02 in		
Height:	65 mm	2.56 in		
Weight max: (incl. contents)	1.94 kg	3.74 lbs		







PNR 6240-87S • First Aid Kit stowed in A320 CAS headrest and hatrack lifeline spacer





EASA AIR-OPS Emergency Medical Kit

- Airbus cross fleet approval
- → ETSO Approval under technical specification TS 40.400/2TS
- → In compliance with EASA AIR-OPS
- → In compliance with FAR 121, subpart X
- In compliance with ICAO, part I
- In case of other needed requirements Innovint can provide a compliance analysis.

Anticipated shelf life is 12 months, which may be reduced down to 9 months due to the availability of medicine production lots.

Latches can be sealed by tie-wraps, plombs or lockers.

The Emergency Medical Kit PNR 6240-22/-62 has been assembled to the statutory requirements of EASA AIR-OPS.

In addition the kit covers FAR 121, subpart X and ICAO, part I as well.

The lower compartment is divided into sections grouping the items by medical subject.

The kit is manufactured of aluminum sheets (coated in gery) and aluminum sections.


The bracket PNR 6240-01-001 allows installation of the kit at any location.

It has been tested in our approved testing facilities for use in all directions (9g).



PNR 6240-22 • Emergency Medical Kit

DIMENSIONS				
Emergency M	led. Kit 🚻	PNR 6240-22		
Emergency M	led. Kit 🧲	PNR 6240-62		
Length: Width: Height:	540 mm 300 mm 115 mm	21.26 in 11.81 in 4.53 in		
Weight max:	6.5 kg	14.32 lbs		



PNR	6240_	.n1_	<i>.</i> 001	Rracket

DIMENSIONS		
Bracket		PNR 6240-01-001
Length:	551 mm	21.69 in
Width:	130 mm	5.12 in
Height:	285 mm	11.22 in
Weight:	1.2 kg	2.64 lbs

ed. 04/15 subject to alterations





FAA Emergency Medical Kit

- Airbus cross fleet approval
- → ETSO Approval under technical specification TS 40.400/2TS
- → In compliance with FAR-121 Subpart X
- In compliance with ICAO, Part I

Anticipated shelf life is 12 months, which may be reduced down to 9 months due to the availability of medicine production lots. Latches can be sealed by tie-wraps, plombs or lockers.

The Emergency Medical Kit PNR 6240-28 has been assembled strictly to the statutory requirements of FAR 121, Subpart X and ICAO, Part I including latest amendment. The emergency medical kit represents the minimum mandatory standard. The upper compartment is equipped with an ampoule rack which holds up to 33 ampoules. The lower compartment is divided into sections grouping the items by medical subject. The kit is manufactured of aluminum sheets (coated in grey) and aluminum sections.


The bracket PNR 6240-01-001

allows installation of the kit at any location. The bracket has been tested in our approved testing facilities for use in all directions (9g).



PNR 6240-28 • Emergency Medical Kit

DIMENSIONS		
Emergency Me	ed. Kit 🔛	PNR 6240-28
Emergency Me	ed.Kit C	PNR 6240-29
Length:	540 mm	21.26 in
Width:	300 mm	11.81 in
Height:	115 mm	4.53 in
Weight max (incl. contents)	5.9 kg	13.00 lbs



PNR	6240-	-01-0i	01 •	Rracket

DIMENSIONS		
Bracket	PNR	6240-01-001
Length:	551 mm	21.69 in
Width: Height:	130 mm 285 mm	5.12 in 11.22 in
Weight:	1,2 kg	2.64 lbs

ed. 04/15 subject to alterations





Attendant Kit

The attendant kit PNR 6240-81 provides for the cabin crew the "all flight" necessary first aid sundry items.

The certified first aid kit does not need to be used, which avoids the removal of the first aid kit from the aircraft and mandatory part inspection.

The soft pocket is made of red Nylon with a zipper. The below listed contents may be extended by some standard medicine, e.g. such as remedy for wounds, nasal decongestant, or loperamide.



PNR 6240-81 /-82 • Attendant Kit

DIMENSIONS		
Attendant Kit		PNR 6240-81
Length: Width: Height:	160 mm 40 mm 110 mm	6.30 in 1.60 in 4.33 in
Weight max:	0,15 kg	0.33 lbs

Attendant Kit	PNR 6240-81

LIST OF CONTENTS

1	EΑ	Adhesive bandage	6 x 10 cm. 8 ea.
2	EΑ	Elastic bandage	8 cm x 4m
2	EΑ	Gauze compress	7,5 x 7,5 cm, 2 ea
1	EΑ	Disposable gloves pair, large	
1	EΑ	Triangular bandage	
1	EΑ	Scissors, bandage	
1	EΑ	Antiseptic alcohol swabs	10 ea

Contents may be extended according to customer requirements, nevertheless space is limited.

DIMENSIONS		
Attendant Kit		PNR 6240-82
Length: Width: Height:	160 mm 40 mm 110 mm	6.30 in 1.60 in 4.53 in
Weight max: (incl. contents)	0,15 kg	0.33 lbs

Same content like PNR 6240-81 with additional medication of paracetamol and gastric anti acid.

ed. 01/15 subject to alterations



ATA CHAPTER 25 60 00

Restraint Kit Standard

- Various possibilities to reliably fix unruly passengers depending on situation
- Includes also means to release fixated persons quickly.

Innovint's Restraint Kit includes all necessary equipment to reliably fixate unruly passengers onboard aircraft.

Various possibilities are given enabling the cabin crew to react adequately: Fixating a passenger on a seat can be done with the heavy duty duct tape or alternatively with the heavy quality chest and leg fixation belts. Passengers lying on their fronts can be bounded either with single use foot shackles and handcuffs or with strong metal handcuffs. To be able to release fixated persons quickly a belt cutter and scissors are included as well.



DIMENSIONS		
Restraint Kit		PNR 7261-07
Length:	250 mm	9.84 in
Width:	80 mm	3.15 in
Height:	110mm	4.33 in
Weight (max):	1.4 kg	3.09 lbs

PNR 7261-07 • Restraint Kit

RESTRAINT KIT

PNR 7261-07 • Restraint Kit

CONTENT

Restraint Kit

PNR 7261-07

1 EA	Heavy duty duct tape
2 EA	Chest/leg fixation belt
2 EA	Single use handcuffs
2 EA	Single use foot shackle
1 EA	Metal handcuffs
1 EA	Belt cutter
1 EA	Scissors





Universal Precaution Kit

Universal Precaution Kits PNR 6240-10X are required on board according to ICAO standard Annex 6 Operation of Aircraft. Innovint adopted the UPKs to its medical emergency equipment product line and offers three different configurations.

The standard UPK is available in a soft robust canvas bag out of the "InnoBag product-line".

This allows an easy stowage throughout the cabin, contributes to a real weight reduction, and fits into an ATLAS standard drawer. Alternatively it can fit in our standard aluminum kit with matching bracket for installation.



PNR 6240-10E • Universal Precaution Kit Standard - soft bag



PNR 6240-10D • Universal Precaution Kit Standard



PNR 6240-10 • Universal Precaution Kit "Long Range", shown with additional lateral hand grip

The PNR 6240-10 is the "Long Range" version and takes, regarding the quantity configuration of the contents, the larger number of crew members into account. Also for the UPK long range version matching installation brackets are available.

On top of the ICAO Standard we offer customized versions, either taking customer requirements into account, or adapt to national regulations regarding the content of the UPK.

Anticipated shelf life is 18 months, which may be reduced down to 15 months due to the availability of medicine production lots.

DIMENSIONS		
UPK Standard		PNR 6240-10 E
Length: Width:	360 mm 240 mm	14.17 in 9.50 in
Height:	95 mm	3.74 in
Weight max: (incl. contents)	2,7 kg	6.00 lbs

DIMENSIONS

UPK Standard	PNR 6240-10D

1	EΑ	Adhesive bandage	6 x 10 cm. 8 ea.
2	EΑ	Elastic bandage	8 cm x 4m
2	EΑ	Gauze compress	7,5 x 7,5 cm, 2 ea
1	EΑ	Disposable gloves pair, large	
1	EΑ	Triangular bandage	
1	EΑ	Scissors, bandage	
1	EΑ	Antiseptic alcohol swabs	10 ea

DIMENSIONS		
UPK Long Rang	ge	PNR 6240-10
Length: Width: Height:	160 mm 40 mm 110 mm	6.30 in 1.60 in 4.33 in
Weight max:	0,15 kg	0.33 lbs

ed. 03/15 subject to alterations

ISSUE DATE: 31.05.10

REV.: N/C

REV. DATE: N/C

1/2

PAGE:

10.

11.

12.

13.

14.

6240-01-362

6240-01-363

6240-01-364

6240-01-365

6240-01-993

LIST OF CONTENTS
INSTRUCTIONS
FOR USE
PNR 6240-10
Universal Precaution Kit



No.	PNR	QTY	Description
1.	6240-01-368	1 x 1 ea.	List of Contents and Instructions for Use in english
2.	6240-01-353	2 x 10 ea.	Face mask (10ea. packing)
3.	6240-01-357	2 x 1 ea.	Eye mask
4.	6240-01-358	2 x 20 ea.	Disposable nitril-gloves (20ea. packing)
5.	6240-01-359	2 x 10 ea.	Protection apron (10ea. packing)
6.	6240-01-355	1 x 100 ea.	Germicidal desinfectant for surface cleaning (100ea. packing)
7.	6240-01-356	4 x 5 ea.	Skin wipes (5ea. packing)
8.	6240-01-354	1 x 150 g	Dry powder (150g)
9.	6240-01-361	1 x 5 ea.	Large absorbent towel, 40cm x 60cm (5ea. packing)

Tie-raps, red (for sealing the box after use)

Biohazard disposable waste bag, 420cm x 610cm (5ea. packing)

Pick up scoop

Scraper with brush

Scraper

1 x 1 ea.

1 x 1 ea.

1 x 1 ea.

1 x 5 ea.

2 x 1 ea.

Abbr.: ea.= each; g= gram

ISSUE DATE: **01. Mar. 2011**

REV.: **N/C**REV. DATE: **N/C**PAGE: **1/2**

LIST OF CONTENTS INSTRUCTIONS FOR USE PNR 6240-10D Universal Precaution Kit



No.	PNR	QTY	Description
1.	6240-01-414	1 x 1 ea.	List of contents and instructions for use; English
2.	6240-01-353	1 x 10 ea.	Face mask (10ea. packing)
3.	6240-01-357	2 x 1 ea.	Eye mask
4.	6240-01-358	1 x 10 pair	Disposable nitril-gloves (10 pair)
5.	6240-01-413	2 x 1 ea.	Protection apron (1ea. packing)
6.	6240-01-356	2 x 5 ea.	Skin wipes (5ea. packing)
7.	6240-01-416	1 x 40 ea.	Germicidal desinfectant for surface cleaning (40ea. packing)
8.	6240-01-362	1 x 1 ea.	Pick up scoop
9.	6240-01-363	1 x 1 ea.	Scraper
10.	6240-01-410	1 x 2 ea.	Large absorbent towel, 40cm x 60cm (2ea. packing)
11.	6240-01-411	1 x 2 ea.	Biohazard disposable waste bag, 420cm x 610cm (2ea. packing)
12.	6240-01-412	2 x 650 g	Dry powder (650g)
13.	6240-01-993	2 x 1 ea.	Tie-raps, red (for sealing the box after use)

Abbr.: ea.= each; g= gram

INNOVINT Aircraft Interior GmbH - Holzmühlenstrasse 84-86 - 22041 Hamburg - Germany - Phone: +49 40 6969860 - Fax: +49 40 6932026 http://www.innovint.de - email: fly@innovint.de

ISSUE DATE: **01. Mar. 2011**REV.: **N/C**

REV. DATE: N/C
PAGE: 1/2

LIST OF CONTENTS INSTRUCTIONS FOR USE PNR 6240-10E Universal Precaution Kit



No.	PNR	QTY	Description
1.	6240-01-415	1 x 1 ea.	List of contents and instructions for use; English
2.	6240-01-353	1 x 10 ea.	Face mask (10ea. packing)
3.	6240-01-357	2 x 1 ea.	Eye mask
4.	6240-01-358	1 x 10 pair	Disposable nitril-gloves (10 pair)
5.	6240-01-413	2 x 1 ea.	Protection apron (1ea. packing)
6.	6240-01-356	2 x 5 ea.	Skin wipes (5ea. packing)
7.	6240-01-416	1 x 40 ea.	Germicidal desinfectant for surface cleaning (40ea. packing)
8.	6240-01-362	1 x 1 ea.	Pick up scoop
9.	6240-01-363	1 x 1 ea.	Scraper
10.	6240-01-410	1 x 2 ea.	Large absorbent towel, 40cm x 60cm (2ea. packing)
11.	6240-01-411	1 x 2 ea.	Biohazard disposable waste bag, 420cm x 610cm (2ea. packing)
12.	6240-01-412	2 x 650 g	Dry powder (650g)
13.	6240-01-993	1 x 1 ea.	Tie-raps, red (for sealing the box after use)

Abbr.: ea.= each; g= gram

Instructions for use
High filtration performance, extremely high absorptive capacity, 95% glass fibre free.
Eye protection glasses to prevent from chemical and mechanical impacts.
Nitril is an artificial rubber material with excellent resistancy against a large number of resolvents and chemicals, as well as oils, fats and acids.
Perfect protection against impureness, made of flexible polyethylene.
Ready-to-use cloths for quick disinfection of surfaces. Medical application: bacteria, fungi and virus to a limited extent.
Reliably acting against vegetative bacteria, mycobacteria, fungi and virus to limited extent.
Highly absorbent inodorous and dust-free granulated material on mineral basis dedicated to absorb liquids.
Cloth to absorb liquids, i.e. blood, secretion or urine, made of cellulose flake with textile fleece sealing. The bottom side consists of a liquid proof polyethylene foil.
Hand shovel according to HACCP (Hazard Analysis and Critical Control Points-conception). Termperature resistant up to 100°C.
Spatula acc. to HACCP (Hazard Analysis and Critical Control Points-
conception). Temperature resistant up to 100°C.